



Model BRW1000 Portable Bariatric Wheelchair Scale Operating Instructions

Caution! This scale weighs approximately 100 pounds. If you think it is too heavy for you to lower or lift safely, ask someone to help you. Damage to the Load Cells can occur if scale is dropped or allowed to fall from the *transport* position. If damage to scale is determined to have been caused by falling or dropping scale, the warranty will be voided!

Lowering Scale from Transport Position to Floor.

1. From transport position, stand close to scale with your feet 8 to 12 inches apart.
2. Firmly grasp handles on scale.
3. While stepping back, lower scale to floor, keeping your back straight and bending your knees until you're in a squatting position.

Never bend at the waist with your legs straight!

DO NOT DROP SCALE!

Lifting Scale from Floor to Transport Position.

1. Stand close to scale with your feet 8 to 12 inches apart.
2. Keeping your back straight, bend your knees until you're in a squatting position.

Never bend at the waist with your legs straight!

3. Firmly grasp handles on scale and lift straight up. Don't twist to either side. Keep scale close to you, not at arm's length.
4. Use your leg muscles as you lift. Keep your back upright and in its natural posture. Lift steadily and smoothly without jerking.

For storage, lay scale flat for maximum protection against damage.



Model BRW1000 Portable Bariatric Wheelchair Scale Operating Instructions

Scale Preparation

1. Install six (6) "C" size batteries or if ordered with scale, plug the small connector end of the AC power adapter cord into the power jack located on the back of the MV1 and then plug the AC power adapter into the proper electrical outlet. **IMPORTANT!** DO NOT connect an AC power adapter to the indicator if alkaline batteries are installed.
2. Install the indicator on the scale bracket and then connect the load cell cable.
3. Place the scale on any hard, level, flat surface or low-cut carpet.
4. The scale is now ready for operation.

NOTE: For more information on installing the batteries (or the AC power adapter), mounting the indicator on the scale bracket, connecting the load cell cable, and setup and operation, refer to the MedVue Model MV1 Medical Weight Analyzer Owner's Manual, 8555-M512-O1.

Basic Operation



WARNING! Only the attendant may move the wheelchair on and off the scale platform. Attendant must assist patient onto and off the scale platform. **NEVER** leave a patient unattended while they are on the scale platform. Failure to maintain control of the patient and/or wheelchair at all times can result in serious injury to you and/or the patient.

To Weigh

1. Press the **ON/OFF** key to turn the MV1 on.
2. Press the **ZERO** key to zero weight display. The ZERO, GROSS and lb or kg annunciator will turn on to show that the scale is ready for use.
3. Assist patient onto scale and read weight display.
4. If a printer is connected to scale, press the **PRINT/ENTER** key to print a ticket.
5. Assist patient off scale.

Zero Weight Display

1. In Gross Weight mode (GROSS annunciator on), press the **ZERO** key.
2. Weight display will return to zero. The ZERO and STABLE **▲▲** annunciators will turn on to show a stable, center-of-zero gross weight condition.

Metric Conversion

Press the **UNITS** key to toggle between pounds and kilograms. The lb or kg annunciator will turn on to show the weighing unit selected. In addition, the height units will change from feet/inches to centimeters. **NOTE:** This is only available if configured in Calibration Setup.



CARDINAL SCALE MFG. CO.
0065-M500-O1 Rev B 203 E. Daugherty, Webb City, MO 64870 USA
06/15 Ph: 417-673-4631 • Fax: 417-673-2153
www.detecto.com

Printed in USA

Technical Support: Ph: 866-254-8261 • tech@cardet.com